

MacroCargo™ Human PBMC-derived Macrophages with IL-11 (Viral System, Lentivirus)

Cat. No.: MTS-1222-YF10

This product is for research use only and is not intended for diagnostic use.

Cell Properties

Product Overview	As a therapeutic tool, macrophage cell has a great capacity for delivering cargos because of their intrinsic characteristics. This product is engineered Human PBMC-derived Macrophages carried with IL-11 by Viral System-Lentivirus. MacroCargo™ products aim to improve the macrophage function and delivery of specific cargos. We also provide custom macrophage delivery systems based on your specific requirements.
Cell Name	PBMC-derived Macrophages
Cell Type	Primary Cell
Cell Species	Human
Cell Background	Macrophages are specialised cells involved in the detection, phagocytosis and destruction of bacteria and other harmful organisms. In addition, they can also present antigens to T cells and initiate inflammation by releasing molecules (known as cytokines) that activate other cells. Furthermore, macrophages are able to cross the biological barriers, penetrate into tumors, and accumulate particularly in hypoxic and poorly vascularized portions of a tumor.

Cargo Properties

Cargo Type	Cytokine
Specific Cargo	IL-11
Cargo Common Name	IL11
Cargo Alternative Names	AGIF; IL-11
Cargo Full Name	Interleukin 11
Introduction	The protein encoded by this gene is a member of the gp130 family of cytokines. These cytokines drive the assembly of multisubunit receptor complexes, all of which contain at least one molecule of the transmembrane signaling receptor IL6ST (gp130). This cytokine is shown to stimulate the T-cell-dependent development of immunoglobulin-producing B cells. It is also found to support the proliferation of hematopoietic

stem cells and megakaryocyte progenitor cells. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

UniprotID	P20809
GeneID	3589
Cargo Delivery System Type	Viral System
Cargo Delivery Approach	Lentivirus

Product Properties

Applications	Improve survive; Promote persistence and activation of endogenous or adoptively transferred NK or T cells
Mycoplasma Testing	Negative
Sterility Testing	Negative
Shipping	Dry ice
Storage	Frozen cells should be stored in a liquid nitrogen tank (-150°C~-190°C) for long term.
Handling Notes	Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance. Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line. Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.
Restriction	Research use only