

MacroCargo™ Human PBMC-derived Macrophages with pDNA encoding IL-5 (Nanoparticle System, Magnetic cationic liposomes)

Cat. No.: MTS-1222-YF25

This product is for research use only and is not intended for diagnostic use.

Cell Properties

Product Overview	As a therapeutic tool, macrophage cell has a great capacity for delivering cargos because of their intrinsic characteristics. This product is engineered Human PBMC-derived Macrophages carried with pDNA encoding IL-5 by Nanoparticle System-Magnetic cationic liposomes. MacroCargo™ products aim to improve the macrophage function and delivery of specific cargos. We also provide custom macrophage delivery systems based on your specific requirements.
Cell Name	PBMC-derived Macrophages
Cell Type	Primary Cell
Cell Species	Human
Cell Background	Macrophages are specialised cells involved in the detection, phagocytosis and destruction of bacteria and other harmful organisms. In addition, they can also present antigens to T cells and initiate inflammation by releasing molecules (known as cytokines) that activate other cells. Furthermore, macrophages are able to cross the biological barriers, penetrate into tumors, and accumulate particularly in hypoxic and poorly vascularized portions of a tumor.

Cargo Properties

Cargo Type	Cytokine
Specific Cargo	pDNA encoding IL-5
Cargo Common Name	IL5
Cargo Alternative Names	EDF; TRF; IL-5
Cargo Full Name	Interleukin 5
Introduction	This gene encodes a cytokine that acts as a growth and differentiation factor for both B cells and eosinophils. The encoded cytokine plays a major role in the regulation of eosinophil formation, maturation, recruitment and survival. The increased production of this cytokine may be related to pathogenesis of eosinophil-dependent inflammatory diseases. This cytokine functions by binding to its receptor, which is a heterod

imer, whose beta subunit is shared with the receptors for interleukine 3 (IL3) and colony stimulating factor 2 (CSF2/GM-CSF). This gene is located on chromosome 5 within a cytokine gene cluster which includes interleukin 4 (IL4), interleukin 13 (IL13), and CSF2 . This gene, IL4, and IL13 may be regulated coordinately by long-range regulatory elements spread over 120 kilobases on chromosome 5q31.

UniprotID	P05113
GeneID	3567
Cargo Delivery System Type	Nanoparticle System
Cargo Delivery Approach	Magnetic cationic liposomes
Nanoparticle Component	DOTAP, DSPC, cholesterol, iron oxide (II, III).

Product Properties

Applications	Improve the delivery of macrophages to tumors and its therapeutic efficacy against inflammatory diseases
Mycoplasma Testing	Negative
Sterility Testing	Negative
Shipping	Dry ice
Storage	Frozen cells should be stored in a liquid nitrogen tank (-150°C~-190°C) for long term.
Handling Notes	Frozen cells should be thawed immediately upon receipt and grown according to handling procedure to ensure cell viability and proper assay performance. Note: Do not freeze the cells upon receipt as it may result in irreversible damage to the cell line. Disclaimer: We cannot guarantee cell viability if the cells are not thawed immediately upon receipt and grown according to handling procedure.
Restriction	Research use only